



### **What must be observed for the shipment of solid wood (unprocessed) to the European Union?**

- Recommended: Phytosanitary certificate. However, Agroisolab is announced as a closed unit for timber as well -> Phytosanitary certificate not always necessary.
- Agroisolab's EU-import number: **DE-NW3 37078** included in the shipment document

### **What is need for test sample analysis?**

- round about >100g; in case of a plywood analysis it is depending on the thickness of the layer (minimum: 1mm). Regularly round about 500g is fine.
- Samples could be dried. Sending in paper envelope or cotton bags recommended.

### **Documentation ?**

- Genus, declared origin
- Proforma invoice including: species, dimension, origin, price